

E-chem E512 pH Meter

Features and benefits

- High precision 3 decimal place benchtop pH meter
- Complete with combination electrode suitable for measurement of low ionic strength solutions.
- User selectable choice of pH buffer calibration sets
- Ability to measure high purity water, suitable for applications in electric power and petrochemical industries.

Technical specifications

Range:

pH:	-2.000 to 19.999 pH
mV:	±1999.9 mV
Temp:	-10 to 110°C

Accuracy:

pH:	±0.002 pH
mV:	±0.03% full scale
Temp:	±0.4°C

Other parameters:

ATC:	0 to 100 °C
Data storage:	600 groups
Comms:	RS 232
Power:	DC9V/300mA
Dimensions:	160x190x70 mm
Weight:	880g

